

Rapid, specific and sensitive visual detection of SARS-CoV-2 (COVID-19) virus RNA using colorimetric duplex loop- mediated isothermal amplification (LAMP)

Disclaimer

This protocol, supported by the [Namur Research Institute for Life Sciences-NARILIS](#) of the [University of Namur](#) (hereafter “we,” “us” or “our”), describes a procedure for SARS-CoV-2 detection. The protocol describes RNA extraction from nasopharyngeal/ oropharyngeal swabs and a duplex LAMP assay.

The provided protocol corresponds to the version valid on June, 2020 and may be subject to changes without notice.

The document is being distributed in good faith and aims at providing a fast and easy to duplicate procedure to scientific community for SARS-CoV-2 detection for diagnosis purposes from clinical specimens. We have taken all reasonable precautions to verify the information contained in this document. However, we make no representation or warranty of any kind, express or implied, regarding the accuracy, adequacy, validity, reliability, availability or completeness of any information therein. The responsibility for the interpretation and use thereof lies with the reader.

Under no circumstance shall we have any liability for any loss or damage of any kind incurred as a result of the use of this document or reliance on any information provided herein.

License

This work is licensed under a [Creative Commons Attribution 4.0 International License](#)



v2020_06_12

Rapid, specific and sensitive visual detection of SARS-CoV-2 (COVID-19) virus RNA using colorimetric duplex loop-mediated isothermal amplification (LAMP)

Keywords : SARS-CoV-2, COVID-19, diagnostic, RNA, LAMP, RT-qPCR, gene N

Axelle Nolmans¹, Damien Coupeau², Benoît Muylkens², Olivier De Backer¹
Contact O. De Backer for enquiries olivier.debacker@unamur.be

¹Namur Research Institute for Life Sciences (NARILIS), Molecular Physiology Research Unit (URPHYM), University of Namur, 61 rue de Bruxelles, 5000 Namur, Belgium.

²Namur Research Institute for Life Sciences (NARILIS), Integrated Veterinary Research Unit (URVI), University of Namur, 61 rue de Bruxelles, 5000 Namur, Belgium.

LAMP is an isothermal nucleic acid amplification technique that allows to detect DNA or RNA molecules of specific sequence such as viral genomes. The reaction is carried at a constant temperature, and does not require expensive thermocyclers used in (q)PCR. It is thus a useful method to perform infectious disease diagnosis at point of care in low- and middle-income countries.

The target sequence is amplified using 2 or more sets of 4 or 6 oligonucleotidic primers each and a polymerase with high strand displacement activity.

The LAMP principle is described here:

<https://www.youtube.com/watch?v=L5zi2P4lggw>

<https://www.youtube.com/watch?v=V4PyySdk3Jo>

LAMP produce amounts of DNA considerably higher than PCR. The amplification product can be detected via different method. The protocol described here is based on a commercial assay (New England Biolabs) and allows visual detection of amplification using pH-sensitive indicator dyes to detect pH change resulting from DNA polymerization.

The protocol was tested on 182 clinical samples which were previously analyzed by conventional RT-qPCR. The duplex LAMP assay detected 100% of the positive samples with $CT \leq 35$ (N=110) in RT-qPCR ; 72 % of samples with $35 < CT < 37$ (N= 18) and a minority of samples with $37 < CT < 40$. The sensitivity was 100% (no false positive among the 27 negative tested samples).

1. Equipment

- Vortex apparatus (only for RNA extraction)
- High speed 1,5ml tube centrifuge (only for RNA extraction)
- Thermostatic oven (preferred) or water bath or dry bath

2. Vessels

- 1.5 tubes for RNA extraction
- 96 well sealable microwell plates and sealing tape (alternatively: 1.5ml or 0.5 ml or 0.2 ml tubes) for LAMP

NB : microwell plates are preferred to tubes because :

- o easy management of a large number of samples.
- o easy record of the visual inspection readout (wells are numbered and a single picture of the plate)
- o possibly more robust results (probably because of reduced volume resulting in less evaporation)

3. Reagents for RNA extraction

- High speed 1,5ml tube centrifuge
- TRI-Reagent from Sigma-Aldrich or TRIzol from Invitrogen or QIAzol from QIAGEN
- Chloroform (Sigma-Aldrich)
- GlycoBlue (ThermoFisher Scientific) (or Glycogen or yeast tRNA or any other RNase free ballast for ethanol co-precipitation to increase visibility of the RNA pellet)
- Isopropanol (Sigma-Aldrich)
- Ethanol 75% (Sigma-Aldrich). In a falcon 50ml tube, add 12,5 ml of RNase free water to 37,5ml of ethanol 96-100% (Sigma-Aldrich).
- RNase free water

4. Reagents for Colorimetric LAMP

- WarmStart® Colorimetric LAMP 2X Master Mix New England Biolabs cat # M1800
- Two mix of 6 primers (synthetic oligonucleotides) each (10 x)

Primers should be HPLC purified (at least FIP and BIP primers)

“GeneN_B” set of primers

GeneN-B-F3	5' ACCGAAGAGCTACCAGACG
GeneN-B-B3	5' TGCAGCATTGTTAGCAGGAT
GeneN-B-FIP	5' TCTGGCCCAGTTCCTAGGTAGTTCGTGGTGGTGACGGTAA
GeneN-B-BIP	5' AGACGGCATCATATGGGTTGCACGGGTGCCAATGTGATCT
GeneN-B-LF	5' CCATCTTGGACTGAGATCTTTCATT
GeneN-B-LB	5' ACTGAGGGAGCCTTGAATACA

“GeneN_ID21” set of primers

F3_ID21	5' GCCAAAAGGCTTCTACGCA
B3_ID21	5' TTGCTCTCAAGCTGGTTCAA
FIP_ID21	5' TCCCCTACTGCTGCCTGGAGGCAGTCAAGCCTCTTCTCG
BIP_ID21	5' TCTCCTGCTAGAATGGCTGGCATCTGTCAAGCAGCAGCAAAG
LF_ID21	5' TGTTGCGACTACGTGATGAGGA
LB_ID21	5' ATGGCGGTGATGCTGCTCT

Prepare the two 10X LAMP Primer Mixtures (one with primers « GeneN-B », the other with primers « ID21 »). Alternatively, a single mixture of all 12 primers can be used.

10X LAMP Primer Mix

FIP	16 µM (1x final concentration in LAMP reaction = 1.6 µM)
BIP	16 µM (1x final concentration in LAMP reaction = 1.6 µM)
F3	2 µM (1x final concentration in LAMP reaction = 0.2 µM)
B3	2 µM (1x final concentration in LAMP reaction = 0.2 µM)
LF	4 µM (1x final concentration in LAMP reaction = 0.4 µM)
LB	4 µM (1x final concentration in LAMP reaction = 0.4 µM)

5. Inactivation of viral infectivity and RNA extraction protocol

Inactivation of viral infectivity and RNA extraction according to Coupeau et al. (https://www.narilis.be/sana/2020_unamur_sars-cov-2_rtqpcr).

Other RNA extraction protocols can be used.

This step and the preparation of LAMP amplification should be performed in a “preLAMP” area to avoid contamination.

RNA is extracted from 100 µl of swab transport medium using TRI Reagent from Sigma-Aldrich or TRIzol from Invitrogen or QIAzol from QIAGEN and dissolved in 30 µl of RNase free water. Store at -20°C (or -80°C if available).

Clinical specimen processing will be performed in a class II biological safety cabinet by experienced researcher wearing appropriate protective equipment (gowns and gloves).

From now, the whole process is carried out using filtered tips.

- Spin the collection tube 200 g, 1min, to make sure the sample is at the bottom of the tube. Use a centrifuge with sealed buckets.
- Clean the centrifuge
- Transfer 100 µl fraction of the clinical sample in a 1,5 ml tube containing 1.0 ml of TRI Reagent.
- Mix gently by inverting the tube. This procedure inactivates sample upon contact with TRI Reagent.
- Incubate 5 min at room temperature. If necessary, inactivated samples can be stored at -20 °C (or -80 °C) for further processing.

- Add 200 µl of chloroform and vortex for 10 sec
- Incubate for 5 min at room temperature
- Centrifuge at 12,000 g for 10 min (preferably at 4 °C)
- Transfer 500 µl of the colorless upper aqueous phase (containing the RNA) in 1.5 ml tube containing 3 µl of Glycoblue (co-precipitant to increase visibility of the RNA pellet)
- Importantly avoid any contact with the ring or the lower organic phase (pink)
- Important note: the organic and aqueous phase can be inverted, i.e. the organic pink phase can be above the clear aqueous phase. In this case, add 100 µl of RNase free water to the sample, vortex and centrifuge at 12,000 g for 10 min preferably at 4 °C)
- Add 500 µl of isopropanol and vortex for 3 sec
- Incubate for 5 min at room temperature
- Centrifuge at 12,000 to 14,000 g for 10 min preferably at 4 °C)
- The RNA pellet forms a blue pellet on the bottom of the tube
- Discard the supernatant
- Add 900 µl of 75 % ethanol
- Gently mix by inverting the tubes (avoid to detach the pellet)
- Centrifuge for 5 min at 12,000 to 14,000 g preferably at 4 °C)
- Aspirate slowly the supernatant with a 1 ml pipette by avoiding contact with the blue pellet.
- Centrifuge briefly again.
- Use a narrow tip to remove residual ethanol
- To dry the pellet, leave the tube open under the chemical hood until complete ethanol evaporation. This step lasts about 5 min. Do not over dry the RNA by letting the sample dry more than 10 min.
- Dissolve the pellet in 30 µl of RNase free water
- Incubate at room temperature until complete dissolution of the blue pellet
- If necessary, RNA can be stored at -80°C for further processing

6. Colorimetric LAMP protocol

Prepare on ice (if available) a mixture (scale up) for the desired number of samples plus positive (positive RNA or synthetic DNA or RNA) and negative controls (negative sample and H₂O).

Per sample:

Colorimetric LAMP 2X Master Mix	12.5 µl
(must be pink when thawed)	
10 X Primer Mix GeneN_B	2.5 µl
10 X Primer Mix GeneN_ID21	2.5 µl
RNase-free H ₂ O	3.5 µl

- Dispense 21 µl in wells (or tubes)
- Put on ice (if available)
- Add 4 µl of RNA sample, mix carefully by up and down pipetting
- Seal the plate with plastic film
- Incubate at 65°C for 35 minutes.

We incubate the plate in a home-made small sandbox placed in an air incubator. Check that the sand is at 65°C.




Sandbox in air incubator.

- Inspect visually the plate.
Positive controls should be yellow, negative controls should be pink. Orange is also to be considered as positive.
- Take a picture of the plate upside down.
See below for an example.
The plate can be frozen and snapped later.
Do not open the plate to avoid contamination.



Electrophoresis readout (optional)

Run 15 µl of the reaction mixture on a 0.7% agarose gel containing ethidium bromide. Positive sample should display amplification products including discrete bands corresponding to the “amplicon” and concatemers of the amplicon (see arrows on the picture).



Gene N sequence showing the location of F3 (sense, left) and B3 (antisense, right) primers as follows:

Gene N-B
ID21

ATGTCTGATAATGGACCCCAAATCAGCGAAATGCACCCCGCATTACGTTTGGTGGACCCTCAGATTC
AACTGGCAGTAACCAGAATGGAGAACGCAGTGGGGCGCGATCAAACAACGTCGGCCCCAAGGTTT
ACCCAATAACTGCGTCTTGGTTCACCGCTCTCACTCAACATGGCAAGGAAGACCTTAAATCCCTC
GAGGACAAGGCGTTCCAATTAACACCAATAGCAGTCCAGATGACCAAATTGGCTACT**ACCGAAGAG**
CTACCAGACGAATTCGTGGTGGTACGGTAAAATGAAAGATCTCAGTCCAAGATGGTATTTCTACTA
CCTAGGAACTGGGCCAGAAGCTGGACTTCCCTATGGTGCTAACAAAGACGGCATCATATGGGTTGCA
ACTGAGGGAGCCTTGAATACACCAAAAAGATCACATTGGCACCCGCA**ATCCTGCTAACAAATGCTGCAA**
TCGTGCTACAACCTTCCCAAGGAACAACATT**GCCAAAAGGCTTCTACGCA**GAAAGGGAGCAGAGGGCG
GCAGTCAAGCCTTCTCGTTCCTCATCACGTAGTCGCAACAGTTCAAGAAATCAACTCCAGGCAGC
AGTAGGGGAACCTTCTCCTGCTAGAATGGCTGGCAATGGCGGTGATGCTGCTCTTGCTTTGCTGCTGC
TTGACAGAT**TTGAACCAGCTTGAGAGCAA**AATGTCTGGTAAAGGCCAACAAACAAGGCCAAACTG
TCACTAAGAAATCTGCTGCTGAGGCTTCTAAGAAGCCTCGGCAAAAACGTAAGTCCACTAAAGCATA
CAATGTAACACAAGCTTTCGGCAGACGTGGTCCAGAACAACCCAAGGAAATTTTGGGGACCAGGA
ACTAATCAGACAAGGAACTGATTACAAACATTGGCCGCAAATTGCACAATTTGCCCCAGCGCTTCA
GCGTCTTCGGAATGTCGCGCATTGGCATGGAAGTCACACCTTCGGGAACGTGGTTGACCTACACAG
GTGCCATCAAATTGGATGACAAAGATCCAAATTTCAAAGATCAAGTCATTTTGCTGAATAAGCATATT
GACGCATACAAAACATTCCCACCAACAGAGCCTAAAAAGGACAAAAAGAAGAAGGCTGATGAAACT
CAAGCCTTACCGCAGAGACAGAAGAAACAGCAAATGTGACTCTTCTCCTGCTGCAGATTTGGATG
ATTTCTCAAACAATTGCAACAATCCATGAGCAGTGCTGACTCAACTCAGGCCTAA